

Head Solutions

2025-2026



P17D730R4



Introduction	4
Head Solution Options	5

i2i	6
-----	---

i2i Overview	7
--------------	---

UniLink™ Ultra	8
----------------	---

Advanced Positioning	10
Highly Adjustable	10
Easy Maintenance	11
How To Order An Ultra	12
Step One	13
Step Two	14
Step Three	15
Step Four	16
Swing Away	17
Alternative Technology	17

Comfort Plus Series	18
---------------------	----

Comfort Plus Configurations	19
CP Series Overview	19
Series Features	20
CP Original	21
Comfort Plus Original size options	21
Fixed Laterals (Tri-Comfort Plus)	22
Tri-Comfort Plus size options	22
Adjustable Lateral Angles	23
Adjustable Comfort Plus size options	23
Modifications	24
Embedded Switches	24
Full Pad Customization	24
CPS Add-Ons	25
Swing Away Receivers	25
Alternative Control Technologies	25
Comparison Chart	26





Combo	27
Combo Overview	28
Simplified Options	29
Swing Away	30
Alternative Technology	30
All Positioning	31
AP Overview	32
AP Pad Sizes	33
Swing Away	34
Alternative Technology	34
Swing Away Receiver	35
SA Overview	36
Alternative Technology Add-Ons	36
Head Solutions Hardware	37
Overview	38
UniLink™	39
Fixed Hardware	39
Flip-Down Hardware	39
HMO Series	40
Basic Hardware (400 Series)	40
500 Series Hardware	40
Ottobock®	41
Accessories	42
Seat Back Mounts (SBM)	42
Adapters	42
OB Push Handle Mount	43
Back Cane Crossbar Mount	44
Ottobock Clamps	44
Protection	45
Warranty Terms	45
Hardware Cover	45
Tone Deflector	45



Head Solutions

How all began

Stealth is a U.S. manufacturer with over 22 years of experience in the rehab market. Launched in 1999, Stealth Products began in response to a need for versatile, high quality head and neck positioning systems, as well as switch site location hardware.

One of our first products introduced to the market, the Stealth Ultra as shown on the right, is a head positioning system that revolutionized the industry.

From our most popular head support, the Comfort Plus, to our most advanced head positioning system, the Unilink™ Ultra, Stealth has a vast experience bringing head positioning solutions to your clients.



Stealth Ultra first generation
launched back in 1999

Head Solution Options



ULTRA



i2i



CP SERIES



AP



CB



Swing Away



Hardware



HEAD SOLUTIONS i2i OVERVIEW

HEAD POSITIONING AT IT'S **FINEST**

The **i2i** is a **head positioning system** that not only **increases proper posture** for the head but for the **upper torso** as well.

The i2i System was developed by Leslie Fitzsimmons PT, ATP who tried all available options for her most challenging students at Lakeview School in New Jersey. **When nothing worked, she designed the i2i System which finally met her client's head positioning needs.**



Small i2i
Low Profile
Short Arms

i2i MAIN FEATURES

UNIQUE SHAPE

The unique, integrated shape of **the i2i provides a foundational connection between the head and neck**, providing support at the **suboccipital shelf** as well as the **occiput**, depending on the profile height. This shape supports the head while limiting neck hyperextension.

ANTERIOR ARMS

Anterior arms promote a midline head position without obstructing shoulder movement. These prevent the client from **hooking** under the head support and also support the **chin prompt**.

Large i2i
Standard Profile
with No Arms and
Ear Cutouts for less
skin friction

Notice that before the i2i the client needed a harness to maintain proper posture

After

Support the suboccipital shelf without contacting the occiput, limiting neck hyperextension.

Before

For the **i2i** complete **detailed** information visit the **official site** or download the **digital publication**





ULTRA



HEAD SOLUTIONS UNILINK™ ULTRA

The Ultra Head Positioning System provides a modular approach to meet the needs of clients with complex presentations.

The occipital pad supports the occiput and is available in a variety of sizes.

The suboccipital pad cups the back of the skull, limiting neck hyperextension, increasing support and stability, and providing lateral support to the neck.

Lateral facial pads can be added to provide lateral support and, in combination with the suboccipital pad, 3 to 4 point control to reduce neck lateral flexion and rotation.




MORE ABOUT
UNILINK ULTRA
HEAD MOUNTS



- Each pad is available in a variety of sizes
- Each pad can be adjusted in multiple planes
- New UniLink mounting hardware is four (4) times holding power stronger than previous hardware at keeping pads in their position
- Foam or gel materials and a variety of upholstery options are available

ULTRA FEATURES HIGHLY ADJUSTABLE



The Ultra's main feature is that it allows double resting point support in the occipital and suboccipital regions.


Both connections support 17mm (0.67in) ball head pads which can be combined to further advance positioning support.

The occipital and suboccipital connections are independent from each other, meaning they can be adjusted to fit more accurately to any asymmetric situation.

ADVANCED POSITIONING

The **Ultra** has been upgraded with the **UniLink** hardware which has **four (4) times** the holding power in each link joint than previous generations. Provides eight main points of adjustments:

- A** Four (4) points of **258°** unrestricted vertical rotation
- B** Two (2) points of **30°** pad angle adjustment and **360°** pad rotation
- C** **6in** (15.24cm) suboccipital horizontal travel length
- D** **360°** of horizontal rotation



MORE ABOUT
UNILINK
HEAD MOUNTS



ULTRA FEATURES EASY MAINTENANCE

REMOVABLE PADS & COVERS

ULTRA STANDARD PADS

The standard Coolcore® covers provide a soft and immersive feel for both pads.. The cover can be removed and replaced. A complete range of custom cover solutions, colors, and styles are available.

Pads are metric

Comes with a standard 5/8in (15.88mm) layer of latex free closed cell rubber foam for support and comfort. It is attached with hook and loop attachment for simple foam replacement and/or maintenance. Gel options are available to reduce friction and increase comfort.

Standard
Coolcore Cover

Latex Free
Closed Cell foam
for support

Hook & Loop layer
for easy removal



MORE ABOUT
COOLCORE



OTHER COVER
OPTIONS

**Pads & Covers
can be customized**

ULTRA ORDERING

HOW TO ORDER AN ULTRA

STEP 2: The **ULTRA** Link

STEP 4: Occipital Pad Options?

STEP 1: How is going to be mounted?


STEP 3: Suboccipital Pad Options?



ULTRA ORDERING STEP ONE

HOW IS IT GOING TO BE MOUNTED?

Want it
Fixed?



Want it
Flip Down?



How much
Height Adjustment?



All hardware comes standard with SBM300-M or SBM306 (15mm) Stealth seat back mount for vertical height adjustment



SBM-DT Dove Tail Mount

Do you need horizontal adjustment?

Only available for 1/2in square staves



Allows for a 2 3/8in (6.03cm) horizontal adjustment

	Item Number	Description
A	UHB2D33200-QCR	Standard 1/2in (12.70mm) square stave fixed hardware with 7in (17.78cm) of height adjustment travel
B	UHB3D33200-QCR	Long 15mm (0.59in) square stave fixed hardware with 12in (30.48cm) of height adjustment travel
C	UHB4D33200-QCR	Long 1/2in (12.70mm) square stave fixed hardware with 12in (30.48cm) of height adjustment travel
D	UHB2D232FD-QCR	Standard 1/2in (12.70mm) square stave flip down hardware with 7 1/2in (19.05cm) of height adjustment travel
E	UHB3D232FD-QCR	15mm (0.59in) square stave flip down hardware with 7 1/2in (19.05cm) of height adjustment travel

All hardware come standard with SBM300M Stealth Seat back mount for vertical height adjustment


ULTRA ORDERING STEP TWO

UHBSU23200 THE ULTRA LINK

This hardware will create a system that utilizes two separate pads to cradle the occipital and suboccipital regions of the head independently. The horizontal bar uses a new clamp style to increase the holding strength and eliminate the potential for damage to the square tube.

Both pad mounts will include a 17mm (0.67in) ball split ring and screws. Along with the mounting hardware the UniLink Ultra hardware has nine (9) adjustment points: Two (2) pad 30° angle and 360° rotation, four (4) vertical 258° angle adjustment, one (1) depth adjustment point for the suboccipital pad, and the height adjustment. All accessible with a single 4mm wrench tool.

If the Dove Tail mount is selected in previous step, then the system will offer ten (10) points of adjustment in one single system.



**Need to add or
replace links?**



ULTRA ORDERING STEP THREE

SUBOCCIPITAL PAD SELECTION



STANDARD SUBOCCIPITAL PADS



	I.N.	Description	Total Width	Maximum Neck Width	Total Height	Contour Depth
A	620M	Large Suboccipital Pad	9½in (24.13cm)	6in (15.24cm)	2¾in (6.99cm)	3½in (8.89cm)
B	600M	Standard Suboccipital Pad	7in (17.78cm)	5in (12.70cm)	3in (7.62cm)	3in (7.62cm)
C	610M	Small Suboccipital Pad	6¼in (15.88cm)	4½in (11.43cm)	2½in (6.35cm)	2¾in (6.99cm)

i2i LOW PROFILE PADS

I.N.	Size	Max Neck Width	Arm Length
I2I250M	Early Intervention	3in (7.62cm)	4 1/2in (11.43cm)
I2I350M	Extra Small	3 1/2in (8.89cm)	4 1/2in (11.43cm)
I2I455M	Small	4in (10.16cm)	6 1/2in (16.51cm)
I2I456M	Small Short Arms	4in (10.16cm)	4 1/2in (11.43cm)
I2I555M	Medium	4 1/2in (11.43cm)	7 1/2in (19.05cm)
I2I556M	Medium Short Arms	4 1/2in (11.43cm)	5 1/2in (13.97cm)
I2I655M	Large	5in (12.70cm)	7 1/2in (19.05cm)
I2I656M	Large Short Arms	5in (12.70cm)	6in (15.24cm)
I2I755M	Extra Large	5 3/4in (14.61cm)	9in (22.86cm)
I2I756M	Extra Large Short Arms	5 3/4in (14.61cm)	6 1/4in (15.88cm)



More **details** about the **i2i**
i2i.stealthproducts.com

ULTRA ORDERING STEP FOUR

OCCIPITAL STANDARD PAD SELECTION

STANDARD OCCIPITAL PADS



	I.N.	Description	Format	Total Width	Total Height	Contour Depth
A	530M	4x2 Occipital Pad	Curved	4in (10.16cm)	2½in (6.35cm)	1/4in (6.35mm)
B	520M	4x3 Occipital Pad	Curved	4in (10.16cm)	3¼in (8.26cm)	1/4in (6.35mm)
C	510M	5x2 Flat Occipital Pad	Flat	5in (12.70cm)	2in (5.08cm)	-
D	500M	7x4 Occipital Pad	Curved	7¼in (18.42cm)	4¼in (10.80cm)	1in (2.54cm)




OTHER OCCIPITAL COMPATIBLE PADS



CP06CCFPM Pediatric Comfort Plus Total Dimensions

6 1/4in (15.88cm) Wide
4 1/4in (10.80cm) Height



CB305M LARGE COMBO Middle Pad Dimensions

5in (12.70cm) Wide
3¼in (8.26cm) Height

CB315M SMALL COMBO Middle Pad Dimensions

5¾in (14.61cm) Wide
4in (10.16cm) Height

ULTRA ADD-ONS SWING AWAY

EXTENDS **ULTRA** POSITIONING & **MORE**



MORE **DETAILS** ABOUT
THE **SWING AWAY RECEIVER**
<https://stlpro.site/stealth-swingaway>



The Ultra is compatible with the **SAR2B** (Low Profile) series. It can be positioned on the occipital, the suboccipital or both. This will allow caretakers or technicians to position any mini-joystick, switches or temporal pads where they are needed most.

Either to enhance positioning or to complement alternative control options, the Swing Away will let a caretaker or technician configure the Ultra to offer solutions to any challenge.


◀ **Shown in the left:** A fiber optic mount with a spot temporal pad swing away arm combo (FO625), red mechanical Egg Switch (ESRED) to change modes, on a low profile left side swing away receiver (SAR2B-L), mounted on the Ultra's occipital region. The Ultra link (UHSU23200) has a large suboccipital pad (620M) and a 7x4 occipital pad (500M), with a standard 1/2in (12.70mm) square tubing FDM hardware (UHB2D232FD-QCR).

ALTERNATIVE TECHNOLOGY

The **UniLink Ultra** can be used as a support for tilt/recline. It can supplement with embedded switches on a tray driving control configuration, or any other configuration, the UniLink Ultra is compatible with our **i-Drive®** Advanced Alternative Control System and **i-Connect™** sensor switches power router.

<https://stlpro.site/i-drive>

<https://stlpro.site/stealth-iconnect>





COMFORT CP PLUS SERIES



HEAD SOLUTIONS

CP SERIES OVERVIEW

The Comfort Plus original is our most popular and affordable contoured head support. The CP Series includes the **Original**, the **Adjustable**, and the **Tri-Comfort Plus**. The generic contour of the **Original CP** encourages a midline head position. The **Fixed** or **Adjustable** angled pads provide more aggressive lateral support, as needed.

The dual layered foam structure reduces pressure when the posterior head is bearing weight, such as during tilt or recline. Wider options can limit **hooking** under the pad.

COMFORT PLUS CONFIGURATIONS

CP


ORIGINAL CONTOURED



Comfort Plus
Original Contour

ACP

ADJUSTABLE COMFORT PLUS



Comfort Plus with
Adjustable Laterals

TCP

TRI-COMFORT PLUS



Comfort Plus with
Fixed Laterals



CP SERIES DETAILS SERIES FEATURES

DUAL LAYERED FOAM RIGID METAL FRAME



Contact layer is a 1½in (3.81cm) open cell contact foam

Base layer is a 5/8in (15.88mm) latex free closed cell foam base

Foam is applied to the base plate with a hook & loop material for easy removal. Making it easy to replace or wash.



COOLCORE® COVER

The CP Series design has been refreshed with a 17mm ball configuration and now comes standard with Coolcore covers.

Coolcore is a research proven performance fabric technology that adapts to its environment and regulates temperature by managing heat and moisture.

The CP Series' covers can be replaced with a wide variety of cover materials to choose from, like our traditional Neoprene cover, Glide Skin and more. Scan OTHER COVER OPTIONS code below to see all available.

Got questions? 1 (800) 965-9229 • 1 (512) 715-9995


A complete range of foams, gels, and memory foams materials are available to add an extra level of personality and comfort. Here are some examples.

Flat Gel


Gel inserts can be added on the inner surfaces for clients who exert strong forces.




SunMate® FRG Foams




Medium Soft




Soft



Extra Extra Soft



Medium



Extra Soft



MORE ABOUT
Coolcore




OTHER COVER
OPTIONS

COMFORT PLUS SERIES CP ORIGINAL

**AN AFFORDABLE
COMFORTABLE PAD**

The **Comfort Plus (CP)** is a single pad headrest made for comfort and ease of use. Its covers and pads are removable and washable. The CP dimensions, covers and pad's foam can be customized to meet the client's requirements.

COMFORT PLUS ORIGINAL SIZE OPTIONS



Customizations

- ✓ Foam Material
- ✓ Cover Material
- ✓ Pad Dimensions

	I.N.	Description	Head Pad Width	Head Pad Height
A	CP06CCFPM	Pediatric Comfort Plus	6¼in (15.88cm)	4¼in (10.80cm)
B	CP08CCFPM	Small Comfort Plus	9¼in (24.29cm)	5¼in (14.61cm)
C	CP10CCFPM	Medium Comfort Plus	11¾in (28.50cm)	5¼in (14.61cm)
D	CP14CCFPM	Large Comfort Plus	14½in (35.96cm)	5¼in (14.61cm)
E	CP19CCFPM	Extra Large Comfort Plus	19¼in (49.53cm)	5¼in (13.02cm)


Part number dimensions refers to base width, *dimensions shown are rounded final dimensions.
This table displays the **true** final dimensions which include base, pad and cover.

COMFORT PLUS SERIES FIXED LATERALS

TRI-COMFORT PLUS 45° FIXED LATERALS

The **Tri-Comfort Plus (TCP)** expands the features of the Comfort Plus with the addition of angled lateral wings and a clinically based choice of length and heights. Increase lateral control with side wings that offer help in keeping the head on the center posterior support.

TRI-COMFORT PLUS SIZE OPTIONS



Customizations

- ✓ Foam Material
- ✓ Cover Material
- ✓ Pad Dimensions

	I.N.	Description	Final Width	Final Height	Inside Length	Wing's Length
A	TCP505M	Pediatric TCP	7¼in (19.69cm)	4½in (11.43cm)	4¾in (12.07cm)	3in (7.62cm)
B	TCP405M	Small TCP	8¼in (22.23cm)	5½in (13.97cm)	6¼in (16.51cm)	3½in (8.89cm)
C	TCP305M	Medium TCP	11½in (29.21cm)	6in (15.24cm)	8½in (21.59cm)	3¾in (9.53cm)
D	TCP205M	Large TCP	14in (35.56cm)	6in (15.24cm)	11in (27.94cm)	5¼in (13.34cm)
E	TCPI05M	Extra Large TCP	18½in (46.99cm)	6in (15.24cm)	15¼in (38.74cm)	7¼in (18.42cm)


*Dimensions shown with pad icons are rounded final dimensions. This table displays the **true** final dimensions which include base, pad and cover.

COMFORT PLUS SERIES ADJUSTABLE

ADJUSTABLE LATERAL ANGLES COMFORT PLUS

The **Adjustable Comfort Plus (ACP)** combines all the features of the Tri-Comfort Plus, adding adjustable lateral angles. Each lateral wing adjusts individually without constraint of any kind.

ADJUSTABLE COMFORT PLUS SIZE OPTIONS



Dimensions shown are with head pad set flat.
Adjust wings as needed.

	I.N.	Description	Final Width	Final Height	Middle Pad Width	Wing's Length
A	ACPI50M	Standard ACP	14in (35.56cm)	5¼in (13.34cm)	5in (12.70cm)	5¼in (13.34cm)
B	ACP550M	Large ACP	19in (48.26cm)	5¼in (13.34cm)	6¾in (17.15cm)	7in (17.78cm)

* Dimensions shown with pad icons are rounded final dimensions. This table displays the **true** final dimensions which include base, pad and cover.

NEED A
LOWER
PROFILE?

Check the
COMBO
Head Pad



1/2in Wing Depth
or Middle Width
Adjustment



INDIVIDUALLY ADJUSTABLE LATERAL ANGLES

Adjust without constraints. Each side
adjust independently.



Customizations

- ✓ Foam Material
- ✓ Cover Material
- ✓ Pad Dimensions

COMFORT PLUS SERIES MODIFICATIONS

**CUSTOMIZATIONS
AVAILABLE FOR CPS**

EMBEDDED SWITCHES



Switches can be embedded into the pad for a clean and discreet look. Either mechanical or proximity switches work with any after market alternative control system or our **Award Winning i-Drive**.

Switches can be used to drive, control tilt/recline or to operate devices like phones and desktops. Switches comes in many formats and shapes like **mechanical** switches and **proximity sensor** switches, for example **Capacitive** or **Fiber Optics**. Some switches cannot be embedded. Contact a representative for more information.

Sip & Puff mechanical switches can be integrated as well with our Swing Away Receiver, Standard or 4 point.




FULL PAD CUSTOMIZATION



The **Original (CP)** and the **Fixed Angles (TCP)** can be customized to fit the customer. From cover, dimensions, and foam to gel additions.

Dimensions provided will be applied to the contour metal frame which will be padded with an approximated **1/4in (6.35mm)** offset around the frame to avoid contact with the cover or client. Contact a Stealth representative or one of our customer service representative members for more information.



Got questions? Call us at one of the phone numbers below

1 (800) 965-9229 or 1 (512) 715-9995

or contact us at information@sw.stealthproducts.com

COMFORT PLUS SERIES CPS ADD-ONS

**CUSTOMIZATIONS
AVAILABLE FOR CPS**

SWING AWAY RECEIVERS



SAR4A-DB






SAR2A-B

The **Comfort Plus Series (CPS)** is compatible with our Swing Away System **Standard (SAR2A)** and **4 point (SAR4A)**. The swing away system helps you place pads, switches, sensor switches, sip & puff, alternative driving methods, or lateral support where it is needed most.



Either to enhance positioning or to complement alternative control options, the Swing Away will let a caretaker or technician configure the Combo to offer solutions to any challenge.

ALTERNATIVE CONTROL TECHNOLOGIES

The **Comfort Plus Series** can be used as a support for tilt or recline. It can supplement with embedded switches on a tray driving control configuration, or any other configuration. It is compatible with our **i-Drive®** Advanced Alternative Control System and **i-Connect™** sensor switches power router.



COMFORT PLUS SERIES COMPARISON CHART



Feature	Comfort Plus	Tri-Comfort Plus	Adjustable Comfort Plus
Format	Curved Frame	Angled Frame (120°)	Tree Part Frame. Center Pad with Adjustable Wings
Positioning Features	Support Only	Aggressive Lateral Support	Custom Lateral Support
Sizes	P S M L XL	P S M L XL	14in (35.56cm) or 19in (48.26cm) Flat
Customizable	Yes	Yes	Yes
Lateral Support	With Swing Away Add-on	Yes	Yes
Swing Away Compatible	Standard (High Profile)		
Foam Structure	Dual Layered		
Tilt / Recline Support	Yes		
Standard Cover	Coolcore		
Cover Options	All Available		
Washable Cover	Yes		
Removable Foam	Yes		

* if further customization needed, contact your Stealth Products' international representatives or our Customer Service Representative.



CB COMBO



HEAD SOLUTIONS COMBO OVERVIEW

The Combo is a triple pad head support that allows movement of the lateral pads from flat (even with the middle pad) to 45°, providing lateral support to the head. The unique contour of the side pads eliminates interference with the ears so that hearing is not impacted and the client's ears do not fold over with head turning. Also, this shape does not interfere with the client's peripheral vision.



**INDIVIDUALLY
ADJUSTABLE
LATERAL PADS**




**LOW PROFILE
CONFIGURATION**



**STANDARD & CUSTOM
LATERAL PADS**

The standard wedge height can be customized upon request to help prevent ear friction against the lateral wing pads.



COMBO OPTIONS SIMPLIFIED OPTIONS




NEED A **LARGER PAD** WITH
ADJUSTABLE LATERALS?
CHECK OUT THE **ACP** AT THE
COMFORT PLUS SERIES



A



B



C

	I.N.	Description	Main Pad Width	Main Pad Height	Wing Pad Width	Wing Pad Height
A	CB315M	Metric Small Combo	5in (12.70cm)	3¼in (8.26cm)	3¼in (8.26cm)	3¼in (8.26cm)
	CB315	Small Combo	5in (12.70cm)	3¼in (8.26cm)	3¼in (8.26cm)	3¼in (8.26cm)
B	CB305M	Metric Large Combo	5½in (14.61cm)	4in (10.16cm)	4¼in (12.07cm)	4in (10.16cm)
	CB305	Large Combo	5½in (14.61cm)	4in (10.16cm)	4¼in (12.07cm)	4in (10.16cm)
C	CB325M	Metric Extra Large Combo	5½in (14.61cm)	6¼in (15.88cm)	7in (17.78cm)	6¼in (15.88cm)
	CB325	Extra Large Combo	5½in (14.61cm)	6¼in (15.88cm)	7in (17.78cm)	6¼in (15.88cm)

Only metric version is compatible with the new UniLink™ head solutions hardware

Adjustable Lateral Hinges



Hinges can extend overall flat pad width **1in (2.54cm)** per side at **1/2in (12.70mm)** increments for a total of **2in (5.08cm)**

Wing Wedge Height



The standard wedge for the small pad is **2in (5.08cm)** and large **2¾in (6.99cm)** but this can be customized at ordering

Removable Pads




5/8in (15.88mm) closed cell latex free foam attached with a **hook & loop** material to the frame for easy removal

COMBO ADD-ONS SWING AWAY

EXTENDS **COMBO**
POSITIONING & **MORE**



MORE **DETAILS** ABOUT
THE **SWING AWAY RECEIVER**
<https://stlpro.site/stealth-swingaway>



SAR2B SERIES
with Arm Height Adjustment

The **metric Combo** is compatible with the **SAR2B** (Low Profile) swing away series while the **original** is compatible with the **SUS9** series. This add-on will allow a caretaker or technician, to position any mini-joystick, switches or temporal pads where they are needed most.

The advantage of the SAR2B series over the SUS9 series is the arm height adjustment feature which has been introduced with the new lighter design.

Either to enhance positioning or to complement alternative control options, the Swing Away will let a caretaker or technician configure the Combo to offer solutions to any challenge.




SUS9 SERIES

ALTERNATIVE TECHNOLOGY

The **Combo** can be used as a support for tilt/recline. It can supplement with embedded switches on a tray driving control configuration, or any other configuration. It is compatible with our **i-Drive®** Advanced Alternative Control System and **i-Connect™** sensor switches power router. The **i-Drive Tri-Array** is based on the Combo.

<https://stlpro.site/i-drive>

<https://stlpro.site/stealth-iconnect>



ALL

POSITIONING



HEAD SOLUTIONS AP OVERVIEW

ALL POSITIONING MALLEABLE HEAD PAD

The **All Positioning (AP)** head support combines occipital and suboccipital support in one head pad. The AP offers a soft padded support for the occipital to rest on during tilt or recline and adjustable arms that can be bent to help provide more lateral suboccipital support. The AP discourages lateral flexion by improving head positioning. Supports swing away system if additional lateral support is required.

DUAL LAYERED FOAM MALLEABLE FRAME



Contact layered is a 1½in (3.81cm)
open cell contact foam

Base layer is a 5/8in (15.88mm)
latex free closed cell foam base

A **malleable** frame that helps to adjust the pad contour to fit client's head contour for better comfort and support

WASHABLE & REMOVABLE

Standard stretchable
Neoprene cover



BASIC POSITIONING



Support to the **occipital** and **suboccipital** with same pad

OCCIPITAL SUPPORT

To **support** during **tilt/recline**


SUBOCCIPITAL SUPPORT

For **positioning** and **comfort**





AP OPTIONS

AP PAD SIZES



Cover can be **customized**
with our variety of options at
www.StealthProducts.com/upholstery



	I.N.	Description	Head Pad Width	Head Pad Height	Suboccipital Concave
A	API50M	Metric Extra Small AP	7in (17.78cm)	6in (15.24cm)	1/2in (12.70mm)
	AP150	Extra Small AP	7in (17.78cm)	6in (15.24cm)	1/2in (12.70mm)
B	AP250M	Metric Small AP	8in (20.32cm)	6¼in (15.88cm)	3/4in (19.05mm)
	AP250	Small AP	8in (20.32cm)	6¼in (15.88cm)	3/4in (19.05mm)
C	AP350M	Metric Medium AP	9in (22.86cm)	6½in (16.51cm)	1in (2.54cm)
	AP350	Medium AP	9in (22.86cm)	6½in (16.51cm)	1in (2.54cm)
D	AP450M	Metric Large AP	10in (25.40cm)	7½in (19.05cm)	1¼in (3.18cm)
	AP450	Large AP	10in (25.40cm)	7½in (19.05cm)	1¼in (3.18cm)
E	AP550M	Metric Extra Large AP	13in (33.02cm)	8in (20.32cm)	1½in (3.81cm)
	AP550	Extra Large AP	13in (33.02cm)	8in (20.32cm)	1½in (3.81cm)

Only metric version is compatible with the new UniLink™ head solutions hardware

AP ADD-ONS SWING AWAY

EXTENDS **HEAD PAD**
POSITIONING & **MORE**



MORE **DETAILS** ABOUT
THE **SWING AWAY RECEIVER**
<https://stlpro.site/stealth-swingaway>



SAR2B SERIES
with Arm Height Adjustment

The metric AP is compatible with the **SAR2B** (Low Profile) swing away series while the original is compatible with the **SUS9** series. The receiver will allow a caretaker or technician to position any mini-joystick, switches or temporal pads where they are needed most.

The advantage of the SAR2B series over the SUS9 series is the arm height adjustment feature which has been introduced with the new lighter design.

Either to enhance positioning or to complement alternative control options, the Swing Away will let a caretaker or technician configure the Combo to offer solutions to any challenge.



SUS9 SERIES



ALTERNATIVE TECHNOLOGY

The **All Positioning** can be used as a support for tilt/recline. It can supplement with embedded switches on a tray driving control configuration, or any other configuration. It is compatible with our **i-Drive®** Advanced Alternative Control System and **i-Connect™** sensor switches power router.

<https://stlpro.site/i-drive>

<https://stlpro.site/stealth-iconnect>

SA SWING AWAY SYSTEM



SWING AWAY RECEIVER SA OVERVIEW

EXTENDING POSITIONING & HEAD ELECTRONICS LOCATION

Our **Swing Away (SA)** receiver supports lateral facial pads, switches for access to external assistive technology devices, or even a chin joystick. The SA has a variety of systems that offer attachment points for support surfaces and switch access.

The SA has **rods** that adjust in angle, rotation, and length. The rods **swing 100° out** of the way for transfers, therapy and when not needed. The rods **swing back** into place **without losing** their **positioning** settings. They also **accommodate** adjustable **temporal pads** which **block lateral flexion** and rotation, and provide a surface for **support**.

LOTS OF
ACCESSORIES
TO CHOOSE
FROM

TEMPORAL PADS
WITH EMBEDDED
SWITCHES

ADDITIONAL RODS
CAN BE ADDED
WITH THE **FIGURE 8
BRACKET** FOR EXTRA
SWITCHES OR JOYSTICK

COVERS OR
FLEX SPIRALS FOR
CLEAN CABLE
MANAGEMENT

360° ROD
ROTATION
120° LATERAL
ROTATION

62° MOUNT
ROTATION

The i-Drive
INTERFACE

SWING AWAY RODS

EASY
SWING AWAY
MECHANISM


ALTERNATIVE TECHNOLOGY ADD-ONS

The **i-Drive®** has integrated the **Swing Away System (SA)** with it's Ultra Pro head control since it's introduction in 2013. This because the SA adds versatility and a high level of customization on-demand. Either implementing the SA as an aftermarket system or integrating it on a fully customized alternative control setting, the SA is a most have for complicated situations.

The SA not only has been used as a Head Solution tool, but has also been implemented on our designs as hand controls positioning mounts as well. The results are a lighter and reliable mounting system that delivers where it is needed the most.

Combining the power of the **i-Drive** and the convenience of the **i-Connect™**, a power router for sensor switches, the SA allows any caretaker or technician to set a control command center that suit the user to achieve the independence this user longs for.

HAVE QUESTIONS? REQUEST THE BROCHURE **HERE**
OR VISIT THE **OFFICIAL WEBSITE**





MHW

MOUNTING HARDWARE



HEAD PAD MOUNTS OVERVIEW

**GREAT HEAD SUPPORT
DESERVES A GOOD FOUNDATION**



17mm Ball Metric Link Style Configuration

The UniLink main focus is simplicity.

One Tool & One System for all configuration.

Fixed
Multi-Axis

Flip-Down
Multi-Axis



17mm Ball & 1in Ball Configuration

HMO Series main focus is affordability.

Strong hardware at the best prices.

Fixed
Basic

Fixed
Offset
Multi-Axis

Flip-Down
Back
Mount



17mm Ball* & 1in Ball Configuration

From a brand that does not need introduction.

Stylish hardware with German engineering.

Fixed
Basic
Multi-Axis

Fixed
Offset
Multi-Axis



*17mm Ball
Adapter



Visit the official website for a full catalog of head solution mounting hardware
<https://stlpro.site/st-hs-hardware>

HEAD PAD MOUNTS UNILINK™

STRENGTH MEETS STYLE
4 TIMES THE HOLDING POWER

FIXED HARDWARE



I.N.	Description	Height Adjustment	Square Stave Area	Back Mount
UHB2D33200	Standard Stave	7in (17.78cm)	1/2in (12.70mm)	SBM300M
UHB3D33200	15mm Long	12in (30.48cm)	15mm	SBM306
UHB4D33200	Standard Long	12in (30.48cm)	1/2in (12.70mm)	SBM300M

The hardware will include split ring, SBM and both metric and SAE fasteners.

FLIP-DOWN HARDWARE



I.N.	Description	Height Adjustment	Square Stave Area	Back Mount
UHB2D232FD	Standard	7 1/2in (19.05cm)	1/2in (12.70mm)	SBM300M
UHB3D232FD	15mm Stave	7 1/2in (19.05cm)	15mm	SBM306

The hardware will include split ring, SBM and both metric and SAE fasteners.

SBM-DT Dove Tail Mount

Do you need horizontal adjustment?

Only available for 1/2in square staves



Allows for a 2 3/8in
(6.03cm) horizontal
adjustment




CHECK ALSO THE NEW
UNILINK ULTRA

The **multi-pad** head positioning
system that revolutionized the
complex rehab industry.



MORE ABOUT
UNILINK
HEAD MOUNTS



HEAD PAD MOUNTS HMO SERIES

AFFORDABLE SOLUTIONS
RELIABLE CONFIGURATIONS

BASIC HARDWARE (400 SERIES)

400 SERIES

I.N.	Height Adjustment	Ball Configuration	Format
HMO400-1	14 $\frac{1}{2}$ in (37.47cm)	1in (2.54cm) Ball	imperial
HMO400-17M	14 $\frac{1}{2}$ in (37.47cm)	17mm	metric

475 / 476 SERIES

I.N.	Height Adjustment	Ball Configuration	Format	Stave Type
HMO475-1	14in (35.56cm)	1in (2.54cm) Ball	imperial	1/2in (1.27cm)
HMO475-17	14in (35.56cm)	17mm	imperial	1/2in (1.27cm)
HMO475-17M	14in (35.56cm)	17mm	metric	1/2in (1.27cm)
HMO476-1	10in (25.40cm)	1in (2.54cm) Ball	imperial	15mm
HMO476-17	10in (25.40cm)	17mm	imperial	15mm

Basic Mounting Hardware (3 Points of Adjustments)

The **400** series is the most basic and affordable option Stealth Products offers. It has 3 points of adjustment: (1) height adjustment, (2) depth adjustment, and (3) pad angle adjustment

Note: Non metric hardware does not include rings.

SURING-SD
17mm Ball
Split Ring



CP-RING-SD
1in Ball
Split Ring

Basic Multi-Axis Mounting Hardware (4 Points of Adjustments)

The **475/476** series is a basic multi-axis hardware system with 4 points of adjustments. (1) It has a vertical angle of 20° and an horizontal 360° radius of adjustment on the triangle. (2) 7in (17.78cm) of depth adjustment. (3) The head pad ball mount allows for a 60° angle of adjustment from any direction and 360° vertical rotation. Finally, (4) height adjustment.

The **475** series hardware includes a **SBM300** back mount. The **476** series includes a **SBM306** back mount. Both metric versions of the hardware series comes with the head pad mounting ring (**SURING-SD**) for 17mm ball configuration.

A **475** series is also available with a **short 10in height adjustment** travel. The **476** series has also a long version available. Call a representative for more information.

500 SERIES HARDWARE

Advanced Multi-Axis Mounting Hardware (7 Points of Adjustments)

I.N.	Height Adjustment	Ball Configuration	Format	Style
HMO500	6in (15.24cm)	universal	metric	standard
HMO501	4in (10.16cm)	universal	metric	pediatric

Mounting packages includes back mount and height adjustment stop rings. Although all rotations are unconstrained, it is recommended not to exceed 45° of rotation.

The **500** series is a low cost multi-axis head support mounting hardware system. It is low profile and light. It has rounded bar for unconstrained and asymmetric adjustments. Comes in standard and pediatric configurations. Another advantage for the 500 series is that is a universal ball mount, which means it support 17mm or 1in (2.54cm) ball head pads.

These series also have 7 points of adjustments: (1) height adjustment, (2) asymmetric vertical stave rotation, (3) 360° vertical bar height rotation, (4) 360° vertical bar angle rotation, (5) depth adjustment, (6) 360° head pad height rotation, and (7) 360° head pad angle rotation.



see it
in action

HEAD PAD MOUNTS OTTOBOCK®

**HEAD MOUNTS FROM A PROVEN
WORLDWIDE RECOGNIZED BRAND**



STANDARD



OFFSET

EXCLUSIVELY DISTRIBUTED BY STEALTH

Proven quality hardware from a leader in the rehabilitation market, OttoBock. 1in (2.54cm) ball configuration, metric. Includes back mount.

Hardware can be adapted to mount 17mm (0.67in) ball head pads with the 17mm adapter ring.

1in Ball Configuration

**HAVE A 17MM BALL PAD?
NO PROBLEM...**



OBSSAR1701
17mm
BALL ADAPTER



OBHDC

OTTOBOCK HEAVY DUTY MOUNTING CLAMP

Only compatible with OttoBock Head Support hardware. The clamp mounts directly into back shell.

I.N.	Description	Height Adjustment	Depth Adjustment	Format
OBHSSMHMSB	Standard (Straight)	11 ¹³ / ₃₂ in (29.01cm)	5 ¹ / ₃₂ in (12.80cm)	metric
OBHSSMHMCB	Offset (Curved)	7 ²³ / ₃₂ in (19.61cm)	5 ¹ / ₃₂ in (12.80cm)	metric



**MORE EXCLUSIVE PRODUCTS
FROM OTTOBOCK BY STEALTH**

<https://stlpro.site/stealth-ob-news>

HEAD PAD MOUNTS ACCESSORIES

HEAD MOUNT ACCESSORIES
EXTEND **HARDWARE** FUNCTIONALITIES

SEAT BACK MOUNTS (SBM)




with 2 $\frac{3}{8}$ in (6.03cm) of
horizontal travel



	I.N.	Description	Format
A	SBM300	1/2in (1.57cm) Square Back Mount	imperial
	SBM300-M	1/2in (1.57cm) Square Back Mount	metric
	SBM306	15mm Square Back Mount	metric
B	SBM-DT	1/2in (1.57cm) Square Dove Tail Mount	metric
C	FDM350	1/2in (1.57cm) Square Flip-Down Back Mount	imperial

ADAPTERS



	I.N.	Description	Format
D	SBM300U	Universal Adapter Plate	imperial
	SBM300-P	Adapter Plate for Permobil® backrests	imperial
	SBM300-KK	Adapter Plate for Kid Kart® backrests	imperial
E	SSAR1701	17mm Adapter Ring for TWB and HMO hardware	metric
	OBSSAR1701	17mm Adapter Ring for Ottobock hardware	metric
F	SBM351	1in (2.54cm) Crossbar Adapter Clamp	imperial
G	HMO-OB	Horizontal adjustment track for SBM or Ottobock hardware	metric

HEAD PAD MOUNTS CANE MOUNTS

WHEELCHAIR MOUNTS FROM OTTOBOCK

OB PUSH HANDLE MOUNT **CLINICAL ASSEMBLY SET**



The Ottobock Push Handle Mount or **Clinical Assembly Set** easily attaches to the horizontal push handles of a wheelchair with a seat width of 15in (38cm) to 24 $\frac{3}{4}$ in (62cm). It is compatible with most headrest mounting hardware. The Clinical Assembly Set fits most wheelchair push handles and its easily installation or removal make it ideal for assessment purposes.

Assembly Set
Main Body

Head Support
Hardware
Track Nut Set

Width
Adjustment
Clamp Bar

CLINICAL ASSEMBLY SET OPTIONS



A



B



C

**EASY TO
REMOVE
NO TOOLS
REQUIRED**



**MORE
FROM
OTTOBOCK
BY STEALTH**

I.N.

Description

A	OBCASMS	Ottobock Clinical Assembly Set for 15in to 16 $\frac{1}{2}$ in (38 - 42cm) wide frames
B	OBCASMM	Ottobock Clinical Assembly Set for 15 $\frac{1}{4}$ in to 21 $\frac{1}{4}$ in (40 - 54cm) wide frames
C	OBCASML	Ottobock Clinical Assembly Set for 20 $\frac{1}{2}$ in to 24 $\frac{1}{4}$ in (52 - 62cm) wide frames

HEAD PAD MOUNTS CANE MOUNTS

WHEELCHAIR MOUNTS FROM OTTOBOCK (Cont.)

The Back Cane Crossbar Mount or **Assembly Set** from Ottobock is a headrest adapter kit which can be added to the vertical cane poles of most wheelchairs. Thanks to the simple mounting system the Assembly Set is removable making it ideal for folding chairs. It is compatible with most headrest mounting hardware.

Fits wheelchair 7 $\frac{1}{8}$ in(18cm) to 24 $\frac{1}{4}$ in(62cm) back frame wide.



BACK CANE CROSSBAR MOUNT ASSEMBLY SET



REQUIRES **TWO (2)**
OTTOBOCK **CLAMPS**
SOLD **SEPARATELY**

	I.N.	Description
A	OBASMS	Ottobock Assembly Set for 7 $\frac{1}{8}$ in to 16 $\frac{1}{2}$ in (18 - 42 cm) wide frames
B	OBASMM	Ottobock Assembly Set for 15 $\frac{3}{4}$ in to 21 $\frac{1}{4}$ in (40 - 54 cm) wide frames
C	OBASML	Ottobock Assembly Set for 20 $\frac{1}{2}$ in to 24 $\frac{1}{4}$ in (52 - 62 cm) wide frames

OTTOBOCK CLAMPS

In most cases these clamps can be attached to the wheelchair without making alterations to the frame. Clamps are size specific.



	I.N.	Description	Profile
	OBSTSCS034	Ottobock Small S.T.S. 3/4in (1.91cm) Clamp	Low
	OBSTSCM078	Ottobock Small S.T.S. 7/8in (2.22cm) Clamp	Low
	OBSTSCL100	Ottobock Small S.T.S. 1in (2.54cm) Clamp	Low
	OBSTSCL078	Ottobock Standard S.T.S 7/8in (2.22cm) Clamp	Standard
	OBC1B	Ottobock Standard S.T.S 1in (2.54cm) Clamp	Standard



**SMALL LOW PROFILE
CLAMPS** SUGGESTED
FOR **OBASMS** CROSSBAR

HEAD PAD MOUNTS PROTECTION

**STEALTH HAS GOT
YOU COVERED!**



MORE INFORMATION ABOUT
STEALTH PRODUCTS' WARRANTY
<https://stlpro.site/returns-n-warranty>

HARDWARE COVER



The **link cover** protects the links from **scratches** and **dust accumulation**. It **also protects** the client by preventing **hair tangling**. Covers can also help with **cable management** on alternative drive control setups.

I.N.	Description
LHW-126	Standard Link Cover for TWB or UniLink hardware
LHW-128	Custom Size Link Hardware Cover


WARRANTY TERMS

**ORDER WITH
CONFIDENCE**



STEALTH PRODUCTS' **LINE OF HARDWARE**
IS **COVERED** BY OUR **5 YEAR**
MANUFACTURING DEFECTS **WARRANTY**

TONE DEFLECTOR



DIFFUSE & ABSORB

If the client is exerting strong forces against the head support, our Tone Deflector can be added to absorb and diffuse these forces, preventing loss of adjustment and breakage.

	I.N.	Description	Format
A	TD100	17mm Ball Tone Deflector	imperial
	TD100M	17mm Ball Tone Deflector	metric
B	TD101	1in (2.54cm) Ball Tone Deflector	imperial

Head Solutions Brochure

Reference Number: P17D730R4

Revision Date: Aug 20, 2025

Customer Satisfaction

Stealth Products is committed to 100% customer satisfaction. Your complete satisfaction is important. Please contact us with feedback or suggested changes that will help improve the quality and usability of our products. You may reach us at:



Stealth Products, LLC

104 John Kelly Drive, Building A, Burnet, Texas 78611

Phone: (512) 715-9995 Toll Free: 1(800) 965-9229

Fax: (512) 715-9954 Toll Free: 1(800) 806-1225

info@stealthproducts.com | www.stealthproducts.com



MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany



Back and **Head Support** sold separately. **Most after-market backs** that can be mounted with a **1in (2.54cm) cane clamp hardware** will fit your new **EPiC Seating**.

Seating setup as shown:

Head Support:
Comfort Plus™

Head Support Hardware:
UniLink™

Back:
ADI® Aluminum Back, also check our TFB™ Backs or ADI Q-Backs

Back Mounting:
Simple Solutions™ FlipLock™

Fits Frame Width: in(cm)

Small 10 to 16 (25.40 to 40.64)

Medium 16 to 21 (40.64 to 53.34)

Large 21 to 26 (53.34 to 66.04)


1in (2.54cm)
Cane Diameter

Quick Adjustment Lever
adjust posterior thoracic support angle without the need of hardware

Gas Struts
come in two resistance levels

50N ideal for clients with limited active movement or posterior pelvic tilt

100N ideal for clients with normal to highly active movement



SCAN ME
FOR MORE
INFORMATION

6in (15.24cm)
posterior pelvic support height*

Pivot Point
EPiC rotates around an imaginary point to the client's acetabulum or the greater trochanter.


* plate height, more options coming soon





Stealth Products

Solutions For Life®



<https://stealthproducts.com>
+1 (800) 965-9229
+1 (512) 715-9995
info@stealthproducts.com

FOLLOW US



104 John Kelly Drive, Building A, Burnet, TX 78654