# JAY FOCUS

The JAY® Focus Point Back is a three-piece contoured backrest featuring a reinforced contoured aluminum shell, multi-adjustable lateral "wings". This back offers a "Lateral Hug" to promote an upright trunk.

Technical Specifications		
Width Range	14" to 21"	
Height	14" to 20"	
Weight	3.5kgs	
Weight Capacity	113kgs	
Avail. Hardware Styles	4-Point	
Angle Adjustability	65°	
Width Adjustability	+1"	
Fits Tubing Sizes	7/8", 1", 1-1/8"	
Headrest Mount Available	Yes	
JAY Your Way Modifications Avail.	No	



# JAY ENCOMPASS

The JAY® Encompass Back is a fully individualised seating solution that offers a wide range of sizes and specifications, including adapted laterals and foam selections. Its three-piece swing-away laterals provide outstanding lateral and anterior stability without inhibiting the ability to transfer.

Technical Specifications	
Width Range	14" to 25"
Height	12" to 28"
Weight	4kgs +
Weight Capacity	113kgs
Avail. Hardware Styles	4-Point
Angle Adjustability	15°
Width Adjustability	N/A
Fits Tubing Sizes	7/8" and 1"
Headrest Mount Available	Yes
JAY Your Way Modifications Avail.	No



# **OPTIONS AND ACCESSORIES**

Please see the product pages for detailed standard availability.

### COVERS

The following cover options are available on select products through the JAY Your Way modification program.

#### SEALED INNER COVER

- Designed to be used underneath outer JAY® cover
- Stretchable Lycra® reduces surface tension
- Polyurethane coating is easy to wipe off and clean
- Aquaguard zipper and anti-wicking seam thread prevent moisture from damaging the foam base

#### JAY STRETCH ZIPPER COVER

- Pleated, 4-way Stretch material reduces surface
- Zipper provides a secure fit and protection for the foam base



#### INCONTINENT RESISTANT SHOWER CAP OR ZIPPER COVER

- Semi-plastic coating can be wiped off and dried to prevent moisture build-up
- Semi-plastic coating prevents moisture from damaging the foam base
- Shower Cap edge makes cover easy to remove
- Zipper provides a secure fit and protection for the foam base



#### X-STATIC® COVER WITH BREATHABLE SPACER FABRIC

- Silver-impregnated material inhibits bacterial growth
- 4-Way stretch reduces surface tension
- Breathable Spacer Fabric promotes airflow



#### DARTEX® COVER WITH AQUAGUARD™ ZIPPER

- Cloth surface prevents sliding
- 4-Way Stretch reduces surface tension
- Dartex® coating and Aquaguard zipper prevent moisture from damaging the foam base



#### REVERSE DARTEX® COVER WITH AQUAGUARD ZIPPER

- Water-resistant Reverse Dartex® coating is easy to wipe off and clean
- 4-Way Stretch reduces surface tension
- Aquaguard zipper prevents moisture from damaging the foam base



# **ACCESSORIES**

#### JAY LUMBAR SUPPORT

- Provides mild support of the lumbar spine
- Helps reduce lower back fatigue and increase sitting tolerance



#### JAY SOLID SEAT

- Improves posture by providing a firm, fully adjustable base of support
- Drop seat option
- Includes mounting hardware





#### **JAY PROTECTOR**

- Provides pressure reduction on a limited basis while the user is sitting or moving about away from the wheelchair
- Fluid-filled cushion helps reduce pressure while foam base assists in stable sitting



## **OPTIONS AND ACCESSORIES**

Please see the product pages for detailed standard availability.

## **BACK HARDWARE**

#### TWO-POINT/SH

- One attachment point on each back cane
- Simple and quick installation
- Angle height & depth adjustable
- Quick release

#### RH / RHF\*

- One attachment point on each back cane
- Lightest hardware solution (quick-release)
- Angle height & depth adjustable
- Quick-release and fixed options available

#### TH\*

- One attachment point on each back cane
- Angle height & depth adjustable
- Lengthened bracket and pin accommodate heavier users or clients with high tone
- Perfect complement for tilt-inspace wheelchairs
- \*Available exclusively on J3 Back.







#### **FOUR-POINT**

- Two attachment points on each back cane
- Withstands greater force and weight by distributing it across a wider area
- Provides rigid mount to wheelchair
- When combined with a Heavy Duty back shell, accommodates increased weight capacities
- Angle height & depth adjustable

#### **STRAP**

- One strap on each side
- Cloth straps click around back posts for easy installation
- Used with the wheelchair's back Upholstery



## LATERAL SUPPORTS

#### **FIXED**

- Provides customized lateral thoracic support
- Width and height adjustable



#### **SWING-AWAY**

- Provides customized lateral thoracic support
- Width and height adjustable
- Swing-away hardware allows for easy transfers out of the wheelchair



#### CRADLE U™

- Provides customized lateral thoracic support for clients with postural or spinal deviations
- Adjusts the existing contour of the JAY J3 back to match the client's unique shape
- Intuitive design attaches directly to the back shell beneath the
- Support plates adjust in angle to match the user's torso



wheelchair Angle adjustability provides the most precise shape to match the user's posture

for easy transfers out of the



