

# Keeping up with changing needs




## Universal Elastic Strap

Temporary positioning for work, play, sports, therapy or fittings with highly stretchable, strong neoprene and extra-fine hook-and-loop that won't snag clothes or trap dirt and lint.

Easy-grasp end tab with a convenient thumb hole for users with limited hand function. Use on the chest, abdomen, lap or legs. The 142cm (56") length can be cut to size without fraying.

- SP110S**      **Small**, 8cm (3") high
- SP110L**      **Large**, 11cm (4-1/2") high
- SP110XL**    **Extra Large**, 15cm (6") high



Attach the thumb loop end of the Universal Elastic Strap to the Extension Strap to for a strong, secure connection.

## Elastic Extension Strap

Combine with the Universal Elastic Strap to increase its length and versatility; for example, to fit around large chairs or carry packages on the lap. Can also be used alone. All extensions are 72.5cm (28") long.

- SP111S**      **Small**, 8cm (3") high
- SP111L**      **Large**, 11cm (4-1/2") high
- SP111XL**    **Extra Large**, 15cm (6") high



Attach and remove effortlessly from the front or side of the wheelchair using the removable plastic seat clip. Alternately, detach the clip and use the hook-and-loop style tabs to attach to an armrest.

## Mobility Bag

The Bodypoint Mobility Bag keeps necessities and valuables close at hand, securely and discreetly. Designed in conjunction with wheelchair users, for wheelchair users, it's big enough to hold the day's needs, but still fits snugly out of the way behind the legs or on the side of the chair.

Sturdy construction and handy features like a two roomy main compartments, a high-contrast interior, mesh front pockets, a key clip, a hanging hook, and easy-grab pulls make this bag unique.

20cm (8") wide, 19cm (7.5") tall, 11cm (4.25") deep

**SP001**

